Determination of Protein Content in Wheat Flour Using Water-cooling-free Fully Automatic Kjeldahl Nitrogen Analyzer

Abstract

Protein content is a key indicator for evaluating the nutritional quality, processing performance and quality control of wheat flour. In accordance with the Kjeldahl method specified in GB 5009.5-2025 National Food Safety Standard — Determination of Protein in Food, this experiment determines the protein reference material in wheat flour with an N810 water-cooling-free fully automatic Kjeldahl nitrogen analyzer, and compares two digestion schemes: powdered catalyst and 15:1 catalyst tablet. The results show that both methods deliver results within the allowable range of reference material with satisfactory accuracy. The catalyst tablet method achieves an RSD of 0.20 %, featuring better precision and simpler operation, which provides data support for the optimization of laboratory detection methods.

Copyright © 2026 Guangzhou Roget Scientific Group Limited All rights reserved.
Telephone: +86-757-83277860 Email: info@rogetscientific.com Website: www.rogetscientific.com
Address: No. 8, Zhixin Road, Dajiao Town, Nansha District, Guangzhou City, Guangdong, China

Introduction

Wheat flour is a staple grain product in the diet of Chinese residents. Protein content serves as a core indicator for the nutritional quality and processing suitability of wheat flour, and is also a critical item for grain quality regulation and pre-delivery inspection in food production.

This experiment is conducted in strict accordance with the first method (Kjeldahl method) of **GB 5009.5-2025** ***National Food Safety Standard — Determination of Protein in Food***. The water-cooling-free fully automatic Kjeldahl nitrogen analyzer is applied to determine the protein reference material in wheat flour. Meanwhile, the detection performance of two digestion schemes, the traditional powdered catalyst system and commercial catalyst tablets, are compared in terms of accuracy, precision and operational convenience, so as to provide data support for the optimization of routine laboratory detection protocols.

Instruments and Reagents

2.1 Instruments

N810 water-cooling-free fully automatic Kjeldahl nitrogen analyzer, NDS-20 Kjeldahl N-digester (20 positions), analytical balance with 0.1 mg readability.

2.2 Reagents

Sulfuric acid (analytical pure), 20 g/L boric acid solution, bromocresol green-methyl red mixed indicator, 400 g/L sodium hydroxide solution, potassium sulfate powder,

Copyright © 2026 Guangzhou Roget Scientific Group Limited All rights reserved.
Telephone: +86-757-83277860 Email: info@rogetscientific.com Website: www.rogetscientific.com
Address: No. 8, Zhixin Road, Dajiao Town, Nansha District, Guangzhou City, Guangdong, China

copper sulfate pentahydrate powder, 15:1 mixed catalyst tablets, 0.1 mol/L hydrochloric acid standard titration solution.

2.3 Sample

Protein reference material (RM) in wheat flour (RM No. BWS0049-2016, Batch No. 2025110706). Certified value: 9.69 g/100g; characteristic value range: 8.91–10.47 g/100g; protein conversion factor: 5.70. (See Annex 1 for the reference material certificate.)

Experimental Procedure

3.1 Sample Preparation

Homogenize the wheat flour protein reference material. Accurately weigh approximately 0.5 g (accurate to 0.1 mg), wrap with weighing paper and transfer into a digestion tube. Prepare reagent blank simultaneously. Two digestion schemes are set as follows:

- **Method 1 (Powdered catalyst)**: Add 6 g of potassium sulfate and 0.4 g of copper sulfate pentahydrate, then slowly add 20 mL of concentrated sulfuric acid along the inner wall of the digestion tube.;
- **Method 2 (Catalyst tablet)**: Add 2 pieces of 15:1 mixed catalyst tablets, then slowly add 10 mL of concentrated sulfuric acid along the inner wall of the digestion tube.

Copyright © 2026 Guangzhou Roget Scientific Group Limited All rights reserved.
Telephone: +86-757-83277860 Email: info@rogetscientific.com Website: www.rogetscientific.com
Address: No. 8, Zhixin Road, Dajiao Town, Nansha District, Guangzhou City, Guangdong, China

Note: To ensure consistent digestion matrix for samples and blanks, a piece of weighing paper is added to the blank digestion tube and digested together with samples.

3.2 Digestion

Place the digestion tubes loaded with samples, catalyst and reagents on the graphite digester. The parameter settings of the digester are listed in Table 1.

Table 1. Digestion Parameter Settings

Stage	Temperature	Time
1	250℃	40min
2	420℃	120min

Copyright © 2026 Guangzhou Roget Scientific Group Limited All rights reserved.
Telephone: +86-757-83277860 Email: info@rogetscientific.com Website: www.rogetscientific.com
Address: No. 8, Zhixin Road, Dajiao Town, Nansha District, Guangzhou City, Guangdong, China

3.3 Instrumental Analysis

After digestion, cool the digestion tubes to room temperature before removal. Check the reagent level of the nitrogen analyzer. Run instrument blanks until three consecutive results are stable, then proceed to sample analysis. The operating parameters of the nitrogen analyzer for Method 1 and Method 2 are shown in Table 2 and Table 3 respectively.

Table 2. Operating Parameters for Method 1

Boric acid volume	*Dilution water volume*	*Alkali solution volume*	*Distillation time*	*Protein conversion factor*	*Titrant concentration*
30mL	30mL	80mL	5min	5.7	0.09998mol/L

Table 3. Operating Parameters for Method 2

Boric acid volume	*Dilution water volume*	*Alkali solution volume*	*Distillation time*	*Protein conversion factor*	*Titrant concentration*
30mL	30mL	40mL	5min	5.7	0.09998mol/L

Copyright © 2026 Guangzhou Roget Scientific Group Limited All rights reserved.
Telephone: +86-757-83277860 Email: info@rogetscientific.com Website: www.rogetscientific.com
Address: No. 8, Zhixin Road, Dajiao Town, Nansha District, Guangzhou City, Guangdong, China

Results and Discussion

4.1 Experimental Results

All samples were processed through digestion, distillation and titration. The test results are summarized in Table 4.

Table 4. Protein Content Test Results

Sample	*Method*	*Sample Weight / g*	*Blank Volume / mL*	*Titration Volume / mL*	*Protein Content (g/100g)*	*Mean (g/100g)*	*Relative Deviation (%)*	*RSD (%)*
Protein quality control sample for wheat flour	Method 1(Powdered catalyst + 20 mL H_2SO_4)	0.5016	0.0248	6.1591	9.76	9.63	-0.62	1.13
		0.5006		6.0304	9.57			
		0.4995		6.0139	9.57			
	Method 2(Catalyst tablet + 10 mL H_2SO_4)	0.5000	0.0373	6.158	9.77	9.79	1.03	0.20
		0.4994		6.1664	9.79			
		0.4997		6.1789	9.81			

4.2 Conclusion

1. Accuracy: The mean protein contents obtained by the two schemes are 9.63 g/100g and 9.79 g/100g respectively, both falling within the characteristic value range of the reference material (8.91–10.47 g/100g). The relative deviations from the certified value (9.69 g/100g) are -0.62 % and 1.03 % respectively, well within the allowable error range of the method. This verifies that both catalytic digestion systems can guarantee the accuracy of detection

Copyright © 2026 Guangzhou Roget Scientific Group Limited All rights reserved.
Telephone: +86-757-83277860 Email: info@rogetscientific.com Website: www.rogetscientific.com
Address: No. 8, Zhixin Road, Dajiao Town, Nansha District, Guangzhou City, Guangdong, China

results.

2. Repeatability: The relative standard deviation (RSD) of Method 1 is 1.13 %, while that of Method 2 is 0.20 %. The catalyst tablet method exhibits superior repeatability, which benefits from the uniform composition of tablet catalysts and high dosing consistency, eliminating operational errors caused by manual powder weighing.

3. Operation and cost: Method 2 only consumes 10 mL of concentrated sulfuric acid, a 50 % reduction compared with Method 1. In addition, catalysts are added directly in tablet form without separate weighing of potassium sulfate and copper sulfate, simplifying sample pretreatment.

References

[1] GB 5009.5-2025, National Food Safety Standard — Determination of Protein in Food[S].

[2] ISO 20483:2013, Cereals and pulses — Determination of the nitrogen content and calculation of the crude protein content — Kjeldahl method[S].

[3] AOAC Official Method 979.09, Protein in Grain, Kjeldahl Method[S].

Copyright © 2026 Guangzhou Roget Scientific Group Limited All rights reserved.
Telephone: +86-757-83277860 Email: info@rogetscientific.com Website: www.rogetscientific.com
Address: No. 8, Zhixin Road, Dajiao Town, Nansha District, Guangzhou City, Guangdong, China

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